

Daily Biologic Calibration Sheet																		
<i>All frequencies must be audible in both ears at 20 dB HL</i>																		
Audiometer Model:										Serial No.:								
Date	Time	Left Ear							Right Ear							Check cords and headphones	Free of Distortion and Unwanted Noise (Dial, Tone button, Cords)	Initials
		250 Hz	500 Hz	1 kHz	2 kHz	4 kHz	6 kHz	8 kHz	250 Hz	500 Hz	1 kHz	2 kHz	4 kHz	6 kHz	8 kHz			

Materials developed by:

Deepa L. Sekhar, M.D., M.Sc.
 Associate Professor of Pediatrics
 Penn State College of Medicine
 Hershey, PA

Beth Czarnecki, Au.D., CCC-A
 Director of Audiology
 Penn State Health
 Milton S. Hershey Medical Center
 Hershey, PA

Thomas R. Zalewski, Ph.D.
 Professor of Audiology
 Bloomsburg University
 Bloomsburg, PA